



**JW Draperies, Inc. Healthcare Services**  
**5625 Furnace Avenue**  
**Elkridge, MD 21075**  
Phone 410-796-9990 Fax 410-796-8912

## **JW *Snap-Mesh* CUBICLE CURTAIN SYSTEM**

The JW Draperies *Snap Mesh* Cubicle Curtain System is standardized to facilitate ease of care and employee safety. The curtain panel and mesh are manufactured separately and then joined by snap tape. To accommodate height variations yet keep the curtains standardized, the height variation is compensated for in the mesh while the fabric curtain panel height remains the same.

### **Mesh Standards**

<b>Fabric: Anti-Microbial</b>	<b>Curtain Height: 64"</b>
<b>8' Ceiling</b>	<b>20" Mesh</b>
<b>8'-6" Ceiling</b>	<b>26" Mesh</b>
<b>9' Ceiling</b>	<b>32" Mesh</b>
<b>10' Ceiling</b>	<b>38" Mesh</b>


Mesh remains in place on the curtain track and the curtain panel only is removed for cleaning. This allows curtain removal without the use of a step ladder and the associated time loss and safety risks.

Standard JW *Snap-Mesh* Curtains are 12' or 18' wide and cover most hospital applications, however curtains are available in any width specified.

Specifying 12' or 18' width curtains is recommended, these sizes will accommodate the majority of hospital curtain tracks and only one panel per track is necessary. Getting the right size curtain in place can be a challenge, the curtain panels are labeled but also, because the curtains are standardized, the width difference is evident by feeling the weight and thickness of the folded curtain.

JW *Snap-Mesh* Curtains come in a standard mesh with integrally woven grommets or are available in a Fabric Border mesh. Both have the Double "Tight Snap" Lead at each end of the curtain, utilizing a tighter lead snap to avoid curtain to mesh disengagement as the curtain is drawn open and closed.

Standardized Snap Count: 6' -18, 9' -26, 12' - 35, 15' -43, 18' -51, 21' -60



**JW Snap-Mesh - Front View - Fabric Border  
Double "Tight Snap" Lead**



JW Snap-Mesh - Rear View - Showing Double "Tight Snap" Leads (Note Gold Band)

## Infection Control Considerations

**Curtains should always be removed for cleaning when soiling is evident, however curtain fabric is available with an Anti-Microbial treatment which may reduce the cleaning frequency requirements (verify with your Infection Control Department).**

**In discussions with Infection Control Managers, and in conjunction with the CDC's position on transmission of disease via fabric, it has been determined that the mesh portion of the JW *Snap-Mesh* Curtain may remain in place because it is above the "touch zone" of the curtain. It is recommended that the Curtain and Mesh panel be removed and cleaned, as a unit, twice annually to maintain a healthy environment (see JW Draperies Healthcare Services Cleaning Program literature).**